

File :Q:\GCMS-06 Minimal\DATA\031120\10-D  
Operator : Organic QC Lab  
Acquired : 11 Mar 2020 12:27 using AcqMethod CHICK\_2019\_VOC-SV\_2S\_2uL.M  
Instrument : GCMS 6  
Sample Name: FRS-054S  
Misc Info in Toluene  
Vial Number: 142



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

