

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-04 Minimal\DATA\032723A\  
... 07- KM.D  
Operator :  
Instrument : GCMS4  
Acquired : 27 Mar 2023 17:05 using AcqMethod VOC-6 SPLIT #H2 17MIN 175DEG.M  
Sample Name: FTOH-004S  
Misc Info : 1H,1H,2H,2H,3H,3H-Perfluorobutan-1-ol in MeOH



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

