

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-04 Minimal\DATA\121423A\
... 57- 223121269 KM.D
Operator :
Instrument : GCMS4
Acquired : 15 Dec 2023 12:06 using AcqMethod VOC-20 SPLIT #H2 15MIN120DEG.M
Sample Name: FTOH-005S
Misc Info : 1H,1H-Perfluorobutan-1-ol in MeOH



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

