

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-04 Minimal\DATA\041323B\  
... 71- KM.D  
Operator :  
Instrument : GCMS4  
Acquired : 14 Apr 2023 19:00 using AcqMethod VOC-20 SPLIT#H2 23MIN.M  
Sample Name: FTOH-015S  
Misc Info : 1H,1H-Perfluorooctan-1-ol in MeOH



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

