

File : Q:\GCMS-04 Minimal\DATA\010524\40- KM.D  
Operator :  
Acquired : 6 Jan 2024 2:29 using AcqMethod VOC-20 SPLIT#H2 28MIN255DEG HOT.M  
Instrument : GCMS4  
Sample Name: FTOH-022S (224011150)  
Misc Info : 1H,1H-Perfluorodecan-1-ol in MeOH  
Vial Number: 11



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

