

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-04 Minimal\DATA\060523\1
... 2- KM.D
Operator :
Instrument : GCMS4
Acquired : 5 Jun 2023 22:01 using AcqMethod VOC-20 SPLIT#H2 23MIN.M
Sample Name: FTOH-026S
Misc Info : 1H,1H-Perfluoroundecan-1-ol in MeOH



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

