

File :K:\QC Data\GCMS-04 Minimal\DATA\010721A\03-RM.D
 Operator : ORGANIC QC DEPT.
 Acquired : 7 Jan 2021 8:54 using AcqMethod VOC-100SPLIT HOT_235 D.M
 Instrument : GCMS4
 Sample Name: GRH-004S/SS
 Misc Info :
 Vial Number: 65

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

