


File : Q:\GCMS-06 Minimal\DATA\091119\02-.D  
Operator : D.C  
Acquired : 11 Sep 19 01:28 pm using AcqMethod CHICK\_2019.M  
Instrument : GCMS 6  
Sample Name : H-112S  
Misc Info : Benz[e]pyrene in Toluene  
Vial Number : 1

 **AccuStandard®**  
Leader in Analytical Reference Standards  
Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
Oven Program: 80c, 17c/min to 340c, 8min  
GC Parameters: Cons. Split, 12psi constant flow  
Split 20:1, 1uL inj.; GC/MS; INJ 270c

