

File :Q:\GCMS-08 Minimal\DATA\20201211B\50- YX.D
Operator : Organic QC Lab
Acquired : 14 Dec 2020 15:00 using AcqMethod Hy_HBDE_S2_2ulinj.M
Instrument : GCMS 8
Sample Name: HBDE-4009S-CN
Misc Info : 5'-Hydroxy-2,3',4,6-tetrabromodiphenyl ether
Vial Number: 115



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

