

File :Q:\GCMS-08 Minimal\DATA\20221219\37- (JP).D
Operator : Organic QC Lab
Acquired : 19 Dec 2022 19:29 using AcqMethod Hy_Pest_24cm_S1.M
Instrument : GCMS 8
Sample Name: HPCB-5010S
Misc Info : 2-Hydroxy-2'34'5'6-Penta-PCB; in iso
Vial Number: 32



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

