

File :Q:\GCMS-08 Minimal\DATA\20191031\70- .D
Operator : KM / JP / ZA
Acquired : 01 Nov 2019 22:14 using AcqMethod KERLA-2017_S50.M
Instrument : GCMS 8
Sample Name: IP-391-07-01
Misc Info : Aromatic Hydrocarbon Test in n-Heptane
Vial Number: 12



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

