

File :Q:\GCMS-06 Minimal\DATA\110519\BSB\18-.D  
Operator : KM  
Acquired : 04 Nov 2019 21:21 using AcqMethod CHICK\_2019\_S100.M  
Instrument : GCMS 6  
Sample Name: ISO-DIS9377-4-2  
Misc Info : Stearyl stearate test solution in Cyclohexane  
Vial Number: 7



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

