

File :L:\QC Data\GCMS-04 Minimal\DATA\081420\05- ZA.D  
Operator : ORGANIC QC DEPT.  
Acquired : 14 Aug 2020 10:22 using AcqMethod VOC-200SPLIT-FAST-11M.M  
Instrument : GCMS4  
Sample Name: M-001R-0.0075X  
Misc Info :  
Vial Number: 1



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

