

File :Q:\GCMS-08 Minimal\DATA\20191127a\J27- .D  
Operator : YX  
Acquired : 28 Nov 19 04:20 am using AcqMethod KERLA-VOC-SV\_50SPLIT.M  
Instrument : GCMS 8  
Sample Name: M-005A-10X  
Misc Info : in DCM  
Vial Number: 59



Column: DB-5, 30m, 0.25 ID, 0.25um  
Oven Program: 40c, 0.5min, 10c/min to 320c, final 2min.  
GC Parameters: Split 50:1; Inlet 250c;

