

File :Q:\GCMS-07 Minimal\DATA\10301901\30- .D  
Operator : KM  
Acquired : 30 Oct 2019 21:10 using AcqMethod KERLA-2017.M  
Instrument : 5977B MSD  
Sample Name: M-1626  
Misc Info Nonylphenol Calibration in MeCl2  
Vial Number: 17



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

