

File :Q:\GCMS-04 Minimal\DATA\112921B\25- YX.D
Operator : Organic QC Lab
Acquired : 30 Nov 2021 19:28 using AcqMethod VOC-50SPLIT DECKER.M
Instrument : GCMS4
Sample Name: M-504.1-CSS
Misc Info : Method 504.1 - Calibration Stock
Vial Number: 61



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

