

File : Q:\GCMS-08 Minimal\DATA\20231122\17- =JP=.D
 Operator : Organic QC Lab
 Acquired : 22 Nov 2023 15:17 using AcqMethod Hy_PAH_15cm_S50.M
 Instrument : GCMS 8
 Sample Name: M-525.3-01
 Misc Info : in acet
 Vial Number: 28



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
 Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
 GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
 Split 2:1, 2µl inj.; GC/MS; Injector 270°C

