

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-06 Minimal\DATA\070822\3
... 1- KM.D
Operator : Organic QC Lab
Instrument : GCMS 6
Acquired : 08 Jul 2022 17:31 using AcqMethod CHICK_2019.M
Sample Name: M-529-SS1
Misc Info : EPA Method 529 Surrogate Standard in MeOH



AccuStandard[®]

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25 μ m
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1 μ l inj.; GC/MS; Injector 270°C

