

File :Q:\GCMS-04 Minimal\DATA\060722\31 RM.D

Operator :

Acquired : 8 Jun 2022 00:21 using AcqMethod VOC-100SPLIT HOT 255 DB624.M

Instrument : GCMS4

Sample Name: M-551.1-IS

Misc Info :

Vial Number: 48



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

