

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-08 Minimal\DATA\20230914
... \19- KM.D
Operator : Organic QC Lab
Instrument : GCMS 8
Acquired : 14 Sep 2023 21:27 using AcqMethod Hy_VOC-SV_15cm_S20.M
Sample Name: M-551.1B
Misc Info : Disinfectant By-Product Standard in Acetone



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

