

File : L:\QC Data\GCMS-04 Minimal\DATA\081220\37 ZA.D
Operator : ORGANIC QC DEPT.
Acquired : 13 Aug 2020 2:31 using AcqMethod VOC-100SPLIT-FAST-11M.M
Instrument : GCMS4
Sample Name: M-551B-4
Misc Info :
Vial Number: 71



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

