

File :Q:\GCMS-08 Minimal\DATA\20190925\19- .D
Operator : ZA / KM / JP
Acquired : 25 Sep 2019 19:06 using AcqMethod KERLA-VOC-SV.M
Instrument : GCMS 8
Sample Name: M-552.1A
Misc Info : HAA Mix in MeOH
Vial Number: 5



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

