

File :Q:\GCMS-08 Minimal\DATA\20190418\17-
Operator : JP/KM
Acquired : 18 Apr 2019 23:12 using AcqMethod KERLA-2017-100uL draw speed.M
Instrument : GCMS 8
Sample Name: M-552.3A-R1 Derivatized
Misc Info : HAA Mix in MtBE Derivatized
Vial Number: 21



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

