

File :L:\QC Data\GCMS-04 Minimal\DATA\062420\167--RM.D
Operator : ORGANIC QC DEPT.
Acquired : 28 Jun 2020 00:12 using AcqMethod VOC-50SPLIT HOT.M
Instrument : GCMS4
Sample Name: M-601-CHG
Misc Info : RT
Vial Number: 3



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

