

File :Q:\GCMS-08 Minimal\DATA\20200123\16- .D
Operator :
Acquired : 23 Jan 2020 16:38 using AcqMethod Hy_Pest_15cm_S20.M
Instrument : GCMS 8
Sample Name:M-622-10
Misc Info : EPN in Hexane
Vial Number: 31



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

