

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-06 Minimal\DATA\052622\3
... 8- KM.D
Operator : Organic QC Lab
Instrument : GCMS 6
Acquired : 27 May 2022 02:17 using AcqMethod VOC-200Split#.M
Sample Name: M-8015B-05
Misc Info : Allyl alcohol in Water



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

