

File : L:\QC Data\GCMS-04 Minimal\DATA\081020\16 ZA.D  
Operator : ORGANIC QC DEPT.  
Acquired : 10 Aug 2020 19:34 using AcqMethod VOC-SCAN 24-50SPLIT-FAST.M  
Instrument : GCMS4  
Sample Name: M-8015B-5031-17  
Misc Info :  
Vial Number: 9



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

