

File :Q:\GCMS-04 Minimal\DATA\080123A\47- RM.D
 Operator :
 Acquired : 2 Aug 2023 9:29 using AcqMethod VOC-100 SPLIT #H2 20MIN175DEG.M
 Instrument : GCMS4
 Sample Name: M-8021B-NAV
 Misc Info : RS
 Vial Number: 83



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm
 Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold
 GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

