

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-08 Minimal\DATA\20231026
 ... \09- KM.D
 Operator : Organic QC Lab
 Instrument : GCMS 8
 Acquired : 26 Oct 2023 13:46 using AcqMethod Hy_Pest_15cm_S20.M
 Sample Name: M-8041
 Misc Info : RCRA Target Phenols Mix in IPA



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
 Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
 GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
 Split 20:1, 1µl inj.; GC/MS; Injector 270°C

