

File :D:\MassHunter\GCMS\1\DATA\080322\06_SAK.D
Operator : Organic QC Lab
Acquired : 03 Aug 2022 09:17 using AcqMethod CHICK_2017_13MIN_S100.M
Instrument : GCMS 6
Sample Name: M-8120-09
Misc Info : 1,2,4,5-Tetrachlorobenzene in Hex
Vial Number: 71



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

