

File : Q:\GCMS-06 Minimal\DATA\071720\19- KM.D  
Operator : Organic QC Lab  
Acquired : 17 Jul 2020 13:08 using AcqMethod CHICK\_2017.M  
Instrument : GCMS 6  
Sample Name: M-8140-20  
Misc Info : Trichloronate in Hexane  
Vial Number: 64



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

