

File :K:\QC Data\GCMS-04 Minimal\DATA\071321\17- RM.D
 Operator :
 Acquired : 13 Jul 2021 19:12 using AcqMethod VOC-100SPLIT DECKER.M
 Instrument : GCMS4 Sample
 Name: M-8240C-R3
 Misc Info :
 Vial Number: 9

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

