

File :Q:\GCMS-08\DATA\20190116\02-  
 Operator : JJP / KM  
 Acquired : 16 Jan 2019 16:22 using AcqMethod KERLA-2017\_S100.M  
 Instrument : GCMS 8  
 Sample Name: M-8270-03  
 Misc Info : SVOCs in MeCl2  
 Vial Number: 77



**AccuStandard®**

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
 Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
 GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
 Split 20:1, 1µl inj.; GC/MS; Injector 270°C

