

File : Q:\GCMS-04 Minimal\DATA\050819\39-.D
Operator : RFXM
Acquired : 9 May 19 4:06 am using AcqMethod VOC-200SPLIT DECKER.M
Instrument : GCMS4
Sample Name: M-GRA-QC-R-10ML
Misc Info :
Vial Number: 31



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

