

File :K:\QC Data\GCMS-04 Minimal\DATA\072921\07- RM.D  
 Operator :  
 Acquired : 29 Jul 2021 13:26 using AcqMethod VOC-SCAN 24 HOT.M  
 Instrument : GCMS4  
 Sample Name: M-GRO-QC-R-IS  
 Misc Info :  
 Vial Number: 61

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

