

File : K:\QC Data\GCMS-04 Minimal\DATA\082420\05--RM.D
Operator : ORGANIC QC DEPT.
Acquired : 24 Aug 2020 21:56 using AcqMethod VOC-100SPLIT DECKER.M
Instrument : GCMS4
Sample Name: M-REF-10
Misc Info :
Vial Number: 41



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

