

File :L:\QC Data\GCMS-04 Minimal\DATA\022820A\44-.D
Operator :
Acquired : 29 Feb 2020 11:34 using AcqMethod VOC-100SPLIT DECKER.M
Instrument : GCMS4
Sample Name: NF-467-CLASS2-B
Misc Info : IN DMSO
Vial Number: 42



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

