

File : Q:\GCMS-08 Minimal\DATA\20200624\41- -DC.D
Operator : Organic QC Lab
Acquired : 25 Jun 2020 00:03 using AcqMethod Hy_Pest_15cm_S10.M
Instrument : GCMS 8
Sample Name: P-031S
Misc Info : Demeton (mixed isomers) in Methanol
Vial Number: 53



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 24c/min to 340c, 2min
GC Parameters: Cons. Split, 5.4psi constant flow
Split 10:1, 1uL inj.; GC/MS; INJ 270c

