

File :D:\MassHunter\GCMS\1\data\11041902\18-.D
Operator : KD&DC
Acquired : 04 Nov 2019 18:54 using AcqMethod KERLA-2017.M
Instrument : 5977B MSD
Sample Name: P-045S-10X
Misc Info in Methanol
Vial Number: 62



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

