

File :D:\MassHunter\GCMS\1\data\20180626\  
Operator : JJP / KM  
Acquired : 26 Jun 2018 22:23 using AcqMethod Pest\_15cm\_S10.M  
Instrument : GCMS 8  
Sample Name: P-080S-10X  
Misc Info : Ronnel 1mg/ml in MeOH  
Vial Number: 23



Column: DB-5, 30m, 0.25 ID, 0.25um  
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

