

File :Q:\GCMS-08\DATA\20190222\04- D  
Operator : JP/KM  
Acquired : 22 Feb 19 10:55 am using AcqMethod KERLA-2017\_S50.M  
Instrument : GCMS 8  
Sample Name: P-125S-H-10X  
Misc Info : Stirophos in Hexane  
Vial Number: 41



**AccuStandard®**

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

