

File :Q:\GCMS-06 Minimal\DATA\021721\26- KM.D
Operator : Organic QC Lab
Acquired : 17 Feb 2021 20:05 using AcqMethod CHICK_2019.M
Instrument : GCMS 6
Sample Name: PLAS-AD-003S
Misc Info : Santoflex IPPD in 80:20 Hexane/A
Vial Number: 47



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

