

File :D:\MassHunter\GCMS\1\data\20200602\11- JP.D
Operator : Organic QC Lab
Acquired : 02 Jun 2020 18:32 using AcqMethod Hy_PAH_15cm_S10.M
Instrument : GCMS 8
Sample Name: R-013S
Misc Info : 3-Nitrofluoranthene in toluene
Vial Number: 57



Column: DB-5, 30m, 0.25ID, 0.25um
Oven: 80c, 2m, 30c/m to 320, final 2min
GC para: Pulse split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj, GCMS, Inj 280c

