

File : Q:\GCMS-08 Minimal\DATA\20200331\15- .D
Operator : Organic QC Lab
Acquired : 31 Mar 2020 18:47 using AcqMethod Hy_PAH_24cm_S100.M
Instrument : GCMS 8
Sample Name: SK-7.5X-100ML
Misc Info : Sulfur in Light distillate kerosene
Vial Number: 79



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 24c/min to 340c, 2min
GC Parameters: Cons. Split, 5.4psi constant flow
Split 100:1 1uL inj.; GC/MS; INJ 270c

