


File :D:\MassHunter\GCMS\1\DATA\010722\14- -KD.D
Operator : Organic QC Lab
Acquired : 07 Jan 22 03:15 pm using AcqMethod CHICK_2017_13MIN_S100.M
Instrument : GCMS 6
Sample Name : UFA-024S
Misc Info : Methyl arachidonate in Heptane
Vial Number: 106

 **AccuStandard**[®]
Leader in Analytical Reference Standards
Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

