

File :Q:\GCMS-08 Minimal\DATA\20231110\33- [JP].D
Operator : Organic QC Lab
Acquired : 10 Nov 2023 21:56 using AcqMethod Hy_PAH_24cm_S200.M
Instrument : GCMS 8
Sample Name: Z-014A
Misc Info : 2mg/ml in dcm
Vial Number: 11



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

