

File :D:\MassHunter\GCMS\1\data\07251801\05  
 Operator : D.C  
 Acquired : 25 Jul 2018 15:44 using AcqMethod CHICK\_2017\_S100.M  
 Instrument : GCMS 6  
 Sample Name: Z-014G-FL  
 Misc Info: Expanded PAH Mix in CH2Cl2:Benzene (50:50)  
 Vial Number: 55

Column: DB-5, 30m, 0.25 ID, 0.25um  
 Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
 GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
 Split 20:1, 1ul inj.; GC/MS; INJ 250c.

