

File :Q:\GCMS-08 Minimal\DATA\20191209\36- .D
Operator : Organic QC Lab
Acquired : 10 Dec 2019 02:19 using AcqMethod KERLA-2017_S100.M
Instrument : GCMS 8
Sample Name: Z-031-13
Misc Info : Simazine in Methanol
Vial Number: 140



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

