

File :D:\MassHunter\GCMS\1\data\09091901\18- .D
 Operator : KD
 Acquired : 09 Sep 2019 15:44 using AcqMethod KERLA-2Split-2uL.M
 Instrument : 5977B MSD
 Sample Name: AE-00050
 Misc Info : in Ethyl acetate
 Vial Number: 77

Column: DB-5, 30m, 0.25 ID, 0.25um
 Oven Program: 100c, 1m, 30c/min to 320c, final 2min.
 GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow
 Split 20:1, 1ul inj.; GC/MS; INJ 250c.

