


File : L:\QC Data\GCMS-06 Minimal\DATA\110519\08-.D
Operator : ZA
Acquired : 05 Nov 2019 14:06 using AcqMethod CHICK_2019_VOC-SV_S100.M
Instrument : GCMS 6
Sample Name: ALR-013S-CN-10X
Misc Info :
Vial Number: 134

 **AccuStandard®**
Leader in Analytical Reference Standards
Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

